



**ASSAY SPECIFICATIONS
FOR**

**FOUNDATIONONE LIQUID CDX
(SPEC-0004, VER. 04)**

**FOUNDATION MEDICINE, INC.
150 SECOND STREET
CAMBRIDGE, MA 02141
USA**

1. Intended Use

FoundationOne Liquid CDx (F1 Liquid CDx) is a next generation sequencing based in vitro diagnostic device that analyzes 324 genes. Substitutions and insertion and deletion alterations (indels) are reported in 311 genes, copy number alterations (CNAs) are reported in 310 genes, and gene rearrangements are reported in 324 genes. The test also detects tumor fraction and the genomic signatures blood tumor mutational burden (bTMB) and microsatellite instability (MSI). F1 Liquid CDx utilizes circulating cell-free DNA (cfDNA) isolated from plasma derived from the anti-coagulated peripheral whole blood of cancer patients. The test is intended to be used as a companion diagnostic to identify patients who may benefit from treatment with targeted therapies in accordance with the approved therapeutic product labeling. Additionally, F1 Liquid CDx is intended to provide tumor mutation profiling to be used by qualified health care professionals in accordance with professional guidelines in oncology for patients with malignant neoplasms.

2. Device Description

The F1 Liquid CDx assay is performed exclusively as a laboratory service using circulating cfDNA isolated from plasma derived from anti-coagulated peripheral whole blood from patients with solid malignant neoplasms. The assay employs a single DNA extraction method to obtain cfDNA from plasma from whole blood. Extracted cfDNA undergoes whole-genome shotgun library construction and hybridization-based capture of 324 cancer-related genes. All coding exons of 309 genes are targeted; select intronic or non-coding regions are targeted in 21 of these genes. Additionally, select intronic or non-coding regions only are targeted in 15 genes (refer to **Table 1** for the complete list of genes targeted by F1 Liquid CDx). Hybrid-capture selected libraries are sequenced with deep coverage using the Illumina NovaSeq® 6000 platform. Sequence data are processed using a custom analysis pipeline designed to accurately detect genomic alterations, including base substitutions, indels, select copy number variants, and select genomic rearrangements. The assay detects substitutions and indels in a total of 311 genes, copy number alterations in 310 genes, and genomic rearrangements in 324 genes. The assay also detects tumor fraction and genomic signatures including MSI and bTMB. A subset of targeted regions in 75 genes is baited for increased sensitivity (**Table 1**).

Table 1: F1 Liquid CDx interrogates 324 genes, including 309 genes with complete exonic (coding) coverage and 15 genes with only select non-coding coverage (indicated with an *); 75 genes (indicated in bold) are captured with increased sensitivity and have complete exonic (coding) coverage unless otherwise noted.

<i>ABL1 [Exons 4-9]</i>	<i>CARD11</i>	<i>DDR1</i>	<i>FGFR3 [Exons 7, 9 (alternative designation exon 10), 14, 18, Intron 17]</i>	<i>KDR</i>	<i>MYD88 [Exon 4]</i>	<i>PPP2R1A</i>	<i>SMARCB1</i>
<i>ACVR1B</i>	<i>CASP8</i>	<i>DDR2 [Exons 5,17,18]</i>	<i>FGFR4</i>	<i>KEAP1</i>	<i>NBN</i>	<i>PPP2R2A</i>	<i>SMO</i>
<i>AKT1 [Exon 3]</i>	<i>CBFB</i>	<i>DIS3</i>	<i>FH</i>	<i>KEL</i>	<i>NF1</i>	<i>PRDM1</i>	<i>SNCAIP</i>
<i>AKT2</i>	<i>CBL</i>	<i>DNMT3A</i>	<i>FLCN</i>	<i>KIT [Exons 8,9,11,12,13,17, Intron 16]</i>	<i>NF2</i>	<i>PRKAR1A</i>	<i>SOCS1</i>

AKT3	CCND1	DOT1L	FLT1	KLHL6	NFE2L2	PRKCI	SOX2
ALK [Exons 20-29, Introns 18,19]	CCND2	EED	FLT3 [Exons 14,15,20]	KMT2A (MLL) [Introns 6, 8-11, Intron 7]	NFKBIA	PTCH1	SOX9
ALOX12B	CCND3	EGFR [Introns 7,15,24-27]	FOXL2	KMT2D (MLL2)	NKX2-1	PTEN	SPEN
AMER1 (FAM123B)	CCNE1	EP300	FUBP1	KRAS	NOTCH1	PTPN11	SPOP
APC	CD22	EPHA3	GABRA6	LTK	NOTCH2 [Intron 26]	PTPRO	SRC
AR	CD70	EPHB1	GATA3	LYN	NOTCH3	QKI	STAG2
ARAF [Exons 4,5,7,11,13,15,16]	CD74* [Introns 6-8]	EPHB4	GATA4	MAF	NPM1 [Exons 4-6,8,10]	RAC1	STAT3
ARFRP1	CD79A	ERBB2	GATA6	MAP2K1 (MEK1) [Exons 2,3]	NRAS [Exons 2,3]	RAD21	STK11
ARID1A	CD79B	ERBB3 [Exons 3,6,7,8,10,12,20,21,23,24,25]	GNA11 [Exons 4,5]	MAP2K2 (MEK2) [Exons 2-4,6,7]	NSD3 (WHSC1L1)	RAD51	SUFU
ASXL1	CD274 (PD-L1)	ERBB4	GNA13	MAP2K4	NT5C2	RAD51B	SYK
ATM	CDC73	ERCC4	GNAQ [Exons 4,5]	MAP3K1	NTRK1 [Exons 14,15, Introns 8-11]	RAD51C	TBX3
ATR	CDH1	ERG	GNAS [Exons 1,8]	MAP3K13	NTRK2 [Intron 12]	RAD51D	TEK
ATRX	CDK12	ERRF1	GRM3	MAPK1	NTRK3 [Exons 16,17]	RAD52	TERC* {ncRNA}
AURKA	CDK4	ESR1 [Exons 4-8]	GSK3B	MCL1	NUTM1* [Intron 1]	RAD54L	TERT* {Promoter}
AURKB	CDK6	ETV4* [Intron 8]	H3F3A	MDM2	P2RY8	RAF1 [Exons 3,4,6,7,10,14,15,17, Introns 4-8]	TET2
AXIN1	CDK8	ETV5* [Introns 6, 7]	HDAC1	MDM4	PALB2	RARA [Intron 2]	TGFBR2
AXL	CDKN1A	ETV6* [Introns 5,6]	HGF	MED12	PARK2	RBI	TIPARP
BAP1	CDKN1B	EWSR1* [Introns 7-13]	HNF1A	MEF2B	PARP1	RBM10	TMPPRSS2* [Introns 1-3]
BARD1	CDKN2A	EZH2 [Exons 4,16,17,18]	HRAS [Exons 2,3]	MEN1	PARP2	REL	TNFAIP3
BCL2	CDKN2B	EZR* [Introns 9-11]	HSD3B1	MERTK	PARP3	RET [Introns 7,8, Exons 11,13-16, Introns 9-11]	TNFRSF14
BCL2L1	CDKN2C	FAM46C	ID3	MET	PAX5	RICTOR	TP53
BCL2L2	CEBPA	FANCA	IDH1 [Exon 4]	MITF	PBRM1	RNF43	TSC1
BCL6	CHEK1	FANCC	IDH2 [Exon 4]	MKNK1	PDCD1 (PD-1)	ROSI [Exons 31,36-38,40, Introns 31-35]	TSC2

<i>BCOR</i>	<i>CHEK2</i>	<i>FANCG</i>	<i>IGF1R</i>	<i>MLH1</i>	<i>PDCD1LG2 (PD-L2)</i>	<i>RPTOR</i>	<i>TYRO3</i>
<i>BCORL1</i>	<i>CIC</i>	<i>FANCL</i>	<i>IKBKE</i>	<i>MPL [Exon 10]</i>	<i>PDGFRA [Exons 12,18, Introns 7, 9, 11]</i>	<i>RSPO2* {Intron 1}</i>	<i>U2AF1</i>
<i>BCR* {Introns 8, 13, 14}</i>	<i>CREBBP</i>	<i>FAS</i>	<i>IKZF1</i>	<i>MRE11A</i>	<i>PDGFRB [Exons 12-21,23]</i>	<i>SDC4* {Intron 2}</i>	<i>VEGFA</i>
<i>BRAF [Exons 11-18, Introns 7-10]</i>	<i>CRKL</i>	<i>FBXW7</i>	<i>INPP4B</i>	<i>MSH2 {Intron 5}</i>	<i>PDK1</i>	<i>SDHA</i>	<i>VHL</i>
<i>BRCA1 {Introns 2, 7, 8, 12, 16, 19, 20}</i>	<i>CSF1R</i>	<i>FGF10</i>	<i>IRF2</i>	<i>MSH3</i>	<i>PIK3C2B</i>	<i>SDHB</i>	<i>WHSC1</i>
<i>BRCA2 {Intron 2}</i>	<i>CSF3R</i>	<i>FGF12</i>	<i>IRF4</i>	<i>MSH6</i>	<i>PIK3C2G</i>	<i>SDHC</i>	<i>WT1</i>
<i>BRD4</i>	<i>CTCF</i>	<i>FGF14</i>	<i>IRS2</i>	<i>MST1R</i>	<i>PIK3CA [Exons 2,3,5-8,10,14,19,21 (Coding Exons 1, 2, 4-7, 9, 13,18,20)]</i>	<i>SDHD</i>	<i>XPO1</i>
<i>BRIP1</i>	<i>CTNNA1</i>	<i>FGF19</i>	<i>JAK1</i>	<i>MTAP</i>	<i>PIK3CB</i>	<i>SETD2</i>	<i>XRCC2</i>
<i>BTG1</i>	<i>CTNNB1 [Exon 3]</i>	<i>FGF23</i>	<i>JAK2 [Exon 14]</i>	<i>MTOR [Exons 19,30,39,40,43-45,47,48,53,56]</i>	<i>PIK3R1</i>	<i>SF3B1</i>	<i>ZNF217</i>
<i>BTG2</i>	<i>CUL3</i>	<i>FGF3</i>	<i>JAK3 [Exons 5,11,12,13,15,6]</i>	<i>MUTYH</i>	<i>PIM1</i>	<i>SGK1</i>	<i>ZNF703</i>
<i>BTK [Exons 2,15]</i>	<i>CUL4A</i>	<i>FGF4</i>	<i>JUN</i>	<i>MYB* {Intron 14}</i>	<i>PMS2</i>	<i>SLC34A2* {Intron 4}</i>	
<i>C11orf30 (EMSY)</i>	<i>CXCR4</i>	<i>FGF6</i>	<i>KDM5A</i>	<i>MYC [Intron 1]</i>	<i>POLD1</i>	<i>SMAD2</i>	
<i>C17orf39 (GID4)</i>	<i>CYP17A1</i>	<i>FGFR1 [Introns 1, 5, Intron 17]</i>	<i>KDM5C</i>	<i>MYCL (MYCL1)</i>	<i>POLE</i>	<i>SMAD4</i>	
<i>CALR</i>	<i>DAXX</i>	<i>FGFR2 [Intron 1, Intron 17]</i>	<i>KDM6A</i>	<i>MYCN</i>	<i>PPARG</i>	<i>SMARCA4</i>	

3. IVDD Product Classification

F1 Liquid CDx is classified in accordance with Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro diagnostic medical devices (IVDD) as a general IVD medical device; it is not listed in Annex II of the IVDD, nor is it a device for self-testing. For general IVD medical devices, the manufacturer self-assesses conformity with the essential requirements and prepares a Declaration of Conformity in accordance with Annex III of the Directive.

F1 Liquid CDx is CE marked according to IVDD 98/79/EC and are registered through Qarad, Foundation Medicine's Authorized Representative in the EU.

4. CE-IVD Part Number

PN-00108

5. Manufacturer

Foundation Medicine, Inc.
150 Second Street
Cambridge, MA 02141
USA



QARAD bvba
Cipalstraat 3, 2440 Geel, Belgium



Foundation Medicine, Inc.
150 Second Street
Cambridge, MA 02141